
Digital Thermometer Instructions – Model DT-105

Congratulations on your purchase of this A&D Medical Digital Fast Read Flex Tip Thermometer. Please read the instructions carefully to ensure accurate temperatures and safe operation. When used properly, this high quality product will give you years of accurate, dependable service.

INDICATIONS FOR USE

The digital thermometer provides a quick and highly accurate reading of an individual's body temperature. The digital thermometer is intended to measure the human body's temperature in regular mode orally, rectally or under the arm, and the device is reusable for clinical or home use on people of all ages.

SAFETY INFORMATION

- ☞ High, prolonged fever requires medical attention (over 100.0°F or 37.8°C) for young children. In case of high readings, please contact your physician.
- ☞ Do not allow children to take their temperature unattended.
- ☞ Read the enclosed instructions carefully to ensure accurate temperature readings.
- ☞ Temperature readings are affected by many factors including exercise, drinking hot or cold beverages prior to taking, and technique.
- ☞ Do not allow children to walk or run during temperature taking. This could present a choking hazard.
- ☞ Oral, underarm, and rectal temperatures vary by approximately 1°F (0.5°C).

SPECIAL PRODUCT FEATURES

- Memory feature: automatic recall of last temperature reading.
- Beeps when peak temperature is reached.
- Auto-Off in approximately 10 minutes.
- Thermometer is waterproof.
- Long-life, replaceable battery (approximately 3,000 temperatures).
- Low Battery Indicator “”.
- Oral, Rectal or Underarm use.

HOW TO CARE FOR THIS PRODUCT

- The thermometer should be cleaned after every use even when probe covers are used.
- Clean the thermometer tip and probe using isopropyl alcohol or soap and water.
- Do not boil the thermometer.
- Store the unit in the protective case while not in use.

- Do not store or keep in direct sunlight, or in high temperatures.
- Do not open the thermometer except to replace the battery.
- Do not drop the thermometer, as this may break the unit.
- Do not attempt to disassemble the unit except to replace the battery.

WHAT IS A “NORMAL” TEMPERATURE?

Temperature readings vary among individuals. A normal temperature refers to each individual's average temperature while they are well. Normal average temperatures range from 97°F/36.1°C to 100°F/37.7°C when taken orally. A rectal temperature is considered your real or “core” body temperature. A rectal temperature is 1°F/0.5°C higher than an oral reading and an underarm temperature is 1°F/0.5°C lower than an oral reading.

HOW TO TAKE TEMPERATURE

1. Select oral, underarm or rectal as your measurement site (see following instructions).
2. Push the ON/OFF button. The thermometer will beep and the display will show all symbols. This is a function check indicating the unit is working properly. If a previous reading is stored in memory, the thermometer will display the last temperature for two seconds. The thermometer is ready for measurement. Begin measurement.
3. When the °F/°C symbol stops flashing and the thermometer will beep for approximately 10 seconds indicating the temperature measurement is complete.
4. Read and record the temperature and time for reference. The reading will automatically be stored in the memory until another reading is taken or the batteries are removed from the thermometer.
5. Push the ON/OFF button to turn the unit off. If you do not turn the unit off, it will shut off automatically in approximately 10 minutes.

NOTE: If measured temperature is less than 90° F (32.2° C), the LCD will display “Lo”. If higher than 109.9° F (42.9° C), “Hi” will be displayed.


ORAL USE

- Do not drink hot or cold fluids, exercise, smoke or perform other activities that raise or lower temperature directly prior to taking reading.
- The mouth should remain closed up to 5 minutes before attempting a reading.
- Place the probe tip well under the tongue in locations as indicated by the “✓” (see diagram). Instruct children not to bite down on the probe.
- The approximate measurement time is 10 seconds. Opening the mouth or improper placement of the probe tip may result in a longer time for a reading.

UNDERARM (AUXILIARY) USE

- From a medical point of view, this method will always give an inaccurate reading and should not be used if a precise measurement is required.
- Wipe underarm with a dry towel.
- Place probe tip in armpit and bring arm down against the body. The measuring sensor should be pointing upward parallel to the arm.
- The approximate measurement time may be 20 seconds or longer depending on the placement of the probe and closure of the arm against the body.

RECTAL USE

- Commonly used for babies, young children or when it is difficult to take an oral or underarm temperature.

- After placing the probe cover on the thermometer, lubricate with a water-soluble jelly such as K-Y Jelly®.
- Gently insert the probe tip no more than 1/2 inch or 1.3 cm into the rectum. Do not force the tip into the rectum if resistance is encountered.
- Wash the thermometer as recommended.
- Approximate measurement time is 8 to 10 seconds.
- Once used rectally, the thermometer should not be used orally for sanitary reasons.

CLEANING INSTRUCTIONS


Wash the area extending about 2” from the tip with warm, soapy water after each use and dry thoroughly.

OR

Wipe the area extending about 2” from the tip with 70% isopropyl alcohol after each use.

CHANGING THE BATTERY

When the “” appears and blinks, the battery needs to be replaced. Make sure you use the same type of battery and proceed as follows: Remove battery compartment lid. Use a non metal probe to remove battery. Insert new equivalent battery type (1.5V or LR41). Ensure battery is correctly positioned with the “+” facing up. Carefully replace the battery compartment lid. Properly dispose of the battery, keeping it from small children and heat.

CHANGING BETWEEN FAHRENHEIT & CELSIUS

1. With the unit off, press and hold the ON/OFF button for approximately 4 seconds to change the current setting.
2. To convert results shown in F° to Celsius (C°) use the following equation:
 $C^{\circ} = 5/9 \times (F^{\circ} - 32)$
3. Here is a helpful chart showing some Fahrenheit and Celsius equivalences:

F°	90.0°	92.0°	94.0°	96.0°	98.0°	100.0°	102.0°	104.0°	106.0°	108.0°
C°	32.2°	33.3°	34.4°	35.6°	36.7°	37.8°	38.9°	40.0°	41.1°	42.2°

SPECIFICATIONS

RANGE	89.6°F - 109.2°F (32.0°C - 42.9°C)
ACCURACY MEASUREMENT RANGE	± 0.2°F for 89.6°F to 107.6°F ± 0.4°F less than 89.6°F ± 0.4°F above 107.6°F ± 0.1°C for 32.0°C to 42.0°C ± 0.2°C less than 32.0°C ± 0.2°C above 42.0°C
STORAGE RANGE	14°F to 131°F (-10°C to 55°C) when stored ≤ 95% RH
DISPLAY	0.1°F (0.1°C)
BATTERY	One 1.5V button size battery
BATTERY LIFE	Approximately 200 hours of continuous operation
DIMENSIONS	4.9 x 0.9 x 0.4” (12.4 x 2.4 x 1.1 cm)
WEIGHT	Approximately 12 grams including battery

LEGAL REQUIREMENTS AND GUIDELINES

This product complies with the European Directive for Medical Devices 93/42/EEC and carries the CE mark. The device also complies with the specifications of the European Standard for clinical thermometers-Part 3:

Performance of compact electrical thermometers (non-predictive and predictive) with maximum device. The CE marking confirms that this is a medical device with a measuring function in the sense of the Medical Devices Act, which has undergone a conformity assessment procedure. A Notified Body confirms that this product fulfills all the appropriate statutory regulations.

COMPLIANCE

Product disposal instructions for electronic devices

The battery in this product complies with the requirements stated in European Directives 2006/66/EEC.

Type B equipment

This device conforms to the European Directive 93/42 EEC for Medical Products.

WARRANTY

A&D Medical warrants to the first retail purchaser that this digital thermometer is free from defects in material and workmanship for two years from date of purchase. This warranty does not cover batteries or damage resulting from improper handling or use. To obtain replacement product, please return the product, freight and insurance prepaid, together with proof of purchase and applicable shipping costs to A&D Medical at the address below. Call the toll-free A&D Medical Health Line at 1-888-726-9966 for shipping costs and further information. At A&D Medical's option, the product will be replaced with a new or refurbished product of equal value. THE WARRANTY AND REMEDY PROVIDED ABOVE ARE EXCLUSIVE AND IN LIEU OF ALL OTHER EXPRESS WARRANTIES. Some states may have laws which amend this warranty.

CONTACT INFORMATION

Our products are designed and manufactured using the latest scientific and technological methods, and offer accurate, easy to use, home monitoring and treatment options. Our full product line includes:

- Digital Blood Pressure Monitors
- Manual Blood Pressure Monitors
- Blood Pressure Cuffs
- Wireless Products
- Personal Health Scales
- Activity Monitors
- Digital Thermometers
- Stethoscopes

This A&D blood pressure product is covered by a two year warranty.

Manufactured for:
A&D Engineering, Inc.
1756 Automation Parkway
San Jose, CA 95131
andmedical.com

Distributed in Canada by:
Auto Control Medical
An A&D Company
6695 Millcreek Drive, Unit 6
Mississauga, Ontario
L5N 5R8 Canada
LifeSourceCanada.com
1-800-461-0991

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