

INSTANT READ Ear Thermometer UT-302

P r e c i s e • I n s t a n t • C o n v e n i e n t



A lightweight and compact thermometer that provides quick and accurate temperature measurements at the push of a button.

The UT-302 Instant Read Ear Thermometer has a unique and attractive design that provides a fast, accurate and easy way to take a temperature through the ear. A temperature reading is provided in just one second! The UT-302 does not require a protective lens cover – its waterproof tip can be easily cleaned with water or rubbing alcohol.

Research indicates that the ear is an ideal site for taking body temperature. The eardrum shares blood vessels with the hypothalamus, the part of the brain that regulates body temperature. Therefore, the ear is an accurate indicator of internal body temperature. Since ear thermometers are designed to receive information from the ear, it is a very safe method for measuring temperature and can be used on children of all ages.

- An accurate reading in just one second (scans the ear 512 times per second)
- Waterproof lens tip that can be easily cleaned with water or rubbing alcohol
- Readings in Fahrenheit or Celsius
- Trilingual instruction guide with illustrations and quick start card in English, French and Spanish

For questions about LifeSource products, call us toll-free at 1-800-461-0991 and a specially trained representative will assist you, or visit us at www.autocontrol.com.

- Fever Indicator
- 10 Memory Recall
- Illuminated Screen
- Battery Included
- One-button Activation
- One-second Reading
- Trilingual Instructions
- 2 Year Warranty

Instant Read UT-302 Ear Thermometer

- Lightweight unit that provides fast and accurate readings.
- Unique shape allows for a gentle and safe method to obtain a temperature reading.
- Perfect for adults and children of all ages.

PART NUMBER	UPC NUMBER	DESCRIPTION
UT-302CN	0-93764-60286-3	Ear Thermometer

Each unit comes with a CR-2032 lithium battery, trilingual instruction manual and quick start card. The UT-302 Ear Thermometer is covered by a two-year limited warranty.

Technical Specifications

Model	UT-302
Dimensions	55.3 x 44.4 x 158.5 mm (2.2" x 1.75" x 6.2")
Weight	65.5 g (2.3 oz.)
Sensor Type	Thermopile Sensor ST118
Measuring range	10°C – 50°C (50°F – 122°F)
Measuring accuracy	± 0.2°C (0.4°F): From 36°C to 39°C (96.8°F - 102.2°F) ± 0.3°C (0.5°F) : Less than 36°C (96.8°F) or greater than 39°C (102.2°F)
Display resolution	0.1°C or °F
Response time	1.0 Second
Memory function	10 memory recall
Thermometer tip	Equipped with a waterproof lens. No lens cover is required.
Automatic shut down	The thermometer will automatically turn off after one minute of non-use.
Power supply	One 3 Volt CR-2032 Lithium battery
Battery life	Approximately 5,000 measurements
Operating environment	16°C – 35°C (60.8°F – 95°F) with relative humidity up to 95%
Storage Environment	-20°C – 50°C (-4°F – 122°F)

Packaging and Shipping Specifications

Selling Unit (Clamshell)	W x D x H	Wt	Individual Item Cube
UT-302	5.9" x 2.0" x 9.1"	0.40 lb	0.06 cu ft
	14.99 cm x 5.08 cm x 23.11 cm	0.18 kg	1.83 cu cm
Unit Case (10 units)	W x D x H	Wt	Case Cube
UT-302	6.2" x 22" x 9.7"	5.0 lb	0.77 cu ft
	15.75 cm x 55.88 cm x 24.64 cm	2.27 kg	23.47 cu cm



Ask your Auto Control Medical Sales Representative about other LifeSource products.

- Blood Pressure Monitors
- Digital Thermometers
- Personal Scales
- Stethoscopes
- Pedometers

For questions about LifeSource products, call us toll-free at 1-800-461-0991 or visit our web site at www.autocontrol.com



6695 Millcreek Drive, Unit 6
Mississauga, ON L5N 5R8
1.800.461.0991 | www.autocontrol.com

