









PULSE OXIMETER FOR EASY AT-HOME MONITORING



KEY FEATURES:

-  Small Size & Lightweight
-  Simple Operation
-  Bar Graph Display
-  Pulse Rate Value Display
-  Low Power Consumption
-  Low Battery Indicator
-  Auto off to conserve battery life
-  2 Year Warranty

Easy to use, reliable Pulse Oximeter

The LifeSource Pulse Oximeter is easy to use with a one button operation. Place the Pulse Oximeter on your finger, press the button, and the results are clearly shown on the LCD. The device is portable making it easy to take measurements on the go and has power saving features to prolong battery life.

The Pulse Oximeter is a non-invasive device intended for checking oxygen saturation of arterial hemoglobin (SpO₂) and the pulse rate of adults at home.

The device can be used for children 4 years and older with a weight greater than 15 kg (33 lbs).

Specifications

| | |
|---------------------------|---|
| Storage Environment | Temperature: -40°C - +60°C Relative humidity: <=95% Atmospheric pressure: 500hPa -1060hPa |
| Operating Environment | Temperature: 10°C - 40°C Relative humidity: <=75% Atmospheric pressure: 700hPa -1060hPa |
| Measurement of SpO2 | Measurement Range: 0 -100% Accuracy: 70-100%, ±2%; Below 70% unspecified |
| Measurement of Pulse Rate | Measurement Range: 30 bpm -250 bpm Accuracy: ±2 bpm |
| Power Requirements | 2 x 1.5V AAA alkaline battery (2.6-3.6V) |
| Power Consumption | Smaller than 25 mA |
| Resolution | SpO2: 1%, Pulse rate: 1 bpm |
| Optical Sensor | Red light (wavelength is 660nm, 6.65mW) Infrared (wavelength is 905nm, 6.75mW) |
| Dimensions | 58.5(L) x 31(W) x 32 (H) mm |
| Weight | ~52g (with the batteries) |

Package Contents

| |
|---------------------|
| Pulse Oximeter |
| Lanyard |
| Carrying Case |
| Batteries (2 x AAA) |

| Wholesaler Code | MSRP |
|------------------|---------|
| McKesson #140209 | \$39.99 |
| LPG #035561 | \$39.99 |

Packaging and Shipping Specifications

| Component | Quantity | Dimensions WxDxH (mm) | Weight (g) | UPC / GTIN |
|----------------|----------|--------------------------|---------------|--------------------|
| Individual Box | 1 | 114.3 x 65.5 x 190 | 127.3 | 0 93764 60394 5 |
| Inner Case | 10 | 460 x 224 x 132 | 1460 | 1 00 93764 60394 2 |
| Outer Case | 100 | 696 x 469 x 496 | 16,200 | 2 00 93764 60394 9 |