









# PULSE OXIMETER FOR EASY AT-HOME MONITORING



## KEY FEATURES:

-  Small Size and Lightweight
-  Auto Rotate
-  Simple Operation
-  Pulse Rate Value Display
-  Low Power Consumption
-  Low Battery Indicator
-  Auto off to conserve battery life
-  2 Year Warranty

## Easy to use, reliable Pulse Oximeter

The A&D Medical Pulse Oximeter is easy to use with a one button operation. Place the Pulse Oximeter on your finger, press the button, and the results are clearly shown in any direction on the LCD. The device is portable making it easy to take measurements on the go and has power saving features to prolong battery life.

The Pulse Oximeter is a non-invasive device intended for checking oxygen saturation of arterial hemoglobin (SpO<sub>2</sub>) and the pulse rate of adults at home.

## Technical Specifications

<b>Storage Environment</b>	Temperature: -40°C - +60°C Relative humidity: <=95% Atmospheric pressure: 500-1060 hPa
<b>Operating Environment</b>	Temperature: 10°C - 40°C Relative humidity: <=75% Atmospheric pressure: 700-106 0hPa
<b>Measurement of SpO2</b>	Measurement Range: 0-100% Accuracy: 70-100%, ±2%; Below 70% unspecified
<b>Measurement of Pulse Rate</b>	Measurement Range: 30-250 bpm Accuracy: ±2 bpm
<b>Power Requirements</b>	2 x 1.5 V AAA alkaline battery (2.6-3.6 V)
<b>Power Consumption</b>	Smaller than 25 mA
<b>Resolution</b>	SpO2: 1%, Pulse rate: 1 bpm
<b>Optical Sensor</b>	Red light (wavelength is 660 nm, 6.65 mW) Infrared (wavelength is 905 nm, 6.75 mW)
<b>Dimensions (LxWxH)</b>	5.9 x 3.3 x 3.2 cm (2.3 x 1.3 x 1.3")
<b>Weight</b>	52 g (1.8 oz) with the batteries

## Packaging and Shipping Specifications

	Quantity	Dimensions (WxLxH)	Weight	Individual Item Cube	UPC / GTIN
<b>Individual Box</b>	1	11 x 4 x 16 cm (4.4 x 1.7 x 6.4")	0.048 kg (1.7 oz)	0.0008 cm <sup>3</sup> (0.03 ft <sup>3</sup> )	0-93764-60415-7
<b>Inner Case</b>	10	34 x 13 x 23 cm (13.3 x 5 x 9.2")	2.5kg (5.5 lbs)	0.01 cm <sup>3</sup> (0.35 ft <sup>3</sup> )	1-00-93764-60415-1
<b>Outer Case</b>	100	68 x 37 x 49 cm (26.6 x 14.4 x 19.2")	18 kg (40 lbs)	0.12 cm <sup>3</sup> (4.26 ft <sup>3</sup> )	2-00-93764-60415-4

## What is In The Box

- 1 Pulse Oximeter
- 1 Instruction Manual
- 1 Carrying Case
- 1 Lanyard
- 2 AAA Batteries

## Wholesaler

	McKesson	LPG	Familiprix	McMahon
<b>Codes</b>	140209	035561	120630	34175501

## Health Canada License - 102228